# GLACE™

Instructions for use

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Bruksanvisning

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Mode d'emploi

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Bruksanvisning

Brugsanvisning

Kasutusjuhised

**Symbols**  
**Symbolit**  
**Symboler**  
**Symbole**  
**Symboles**  
**Symbolen**  
**Simboli**  
**Symboler**  
**Symboler**  
**Sümbolid**

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
## ENGLISH

### **Description**

The GLACE™ Custom made craniomaxillofacial Implant is a STERILE, single use plate intended for use in craniomaxillofacial surgery. The plate is bacteriostatic and promotes bone ingrowth due to the bioactive glass incorporated into the plate. The plate is surgically implanted and fixed in place with titanium screws. All GLACE™ Custom Made CMF Implant designs have been approved by the implanting surgeon prior to manufacture of the device.

Typical usage situations are: reconstruction of cranial bone defects and damage (e.g. resulting from trauma, neoplasm resection, infected craniotomy flaps or neurosurgical operations), reconstruction of cranial and maxillofacial bone regions during surgery, including reconstruction of orbital floor.

### **Indications**

The intended use for the GLACE™ Custom Made CMF Implant is augmentation and reconstruction of the craniomaxillofacial skeleton. It is indicated for non-load bearing applications for all patient age groups, and for use with an intact dura, with or without duraplasty.

### **Contraindications**

- Active infection.
- Foreign body sensitivity. Where material sensitivity is suspected, testing is to be completed prior to implantation.
- Patients with limited blood supply, insufficient quantity or quality of bone, insufficient soft tissue quantity or quality, or latent infection.
- Patients with mental or neurological conditions who are unwilling or incapable of following postoperative care instructions.

### **Instructions for clinical use**

- Read these instructions for use in their entirety before beginning surgery.
- Confirm the GLACE™ Custom Made CMF Implant product label for: patient ID, expiration date.
- Confirm the GLACE™ Custom Made CMF Implant intact sterility barrier before beginning surgery or opening the GLACE™ Custom Made CMF Implant packaging.
- The GLACE™ Custom Made CMF Implant is manufactured to fit the patient's defect site. The implant is designed to exceed the region of the bone defect by

3–20 mm. In order to allow the implant to be fixed over the bone, the implant edge is minimum 6 mm. An edge of less than 6 mm does not have pre-made screw holes and the manufacturer does not recommend drilling additional holes to these edges.

- Reconstruction can be made with several implants which overlap each other, but all of them must be fixed to the remaining bone.
- The fixation screws are installed through the holes at the edges of the implant. We recommend flat head self-drilling titanium screws for the attachment, but other screws may be used. The number of screws should be determined by the operating surgeon based on the size and shape of the implant.
- The GLACE™ Custom Made CMF Implant allows minor bending by tightening the screws. Please take care that the bending does not break the composite structure.
- The implant's structure is porous and intended to absorb blood during operation. To optimise blood absorption, the implant must be implanted dry.
- The manufacturer recommends that the top of the skull next to defect margins be refreshed. The dura mater should be protected during refreshing. A dura mater tack-up suture can be implemented using the tack-up suture holes in the implant or holes separately commissioned by the surgeon in conjunction with the order. The cranial flap is closed and sutured carefully to cover the implant.

### **Possible adverse effects and complications**

1. Poor bone formation, osteoporosis, osteolysis, osteomyelitis, inhibited revascularisation, or infection can cause the device to loosen or fracture.
2. Non-union or delayed union, which may lead to the displacement, loosening or breakage of the implant.
3. Sensitivity to device materials, or allergic reaction to a foreign body.
4. Pain, discomfort, abnormal sensation, or palpability due to the presence of the device.
5. Increased fibrous tissue response around the implant.
6. Inadequate healing.
7. Hematoma and seroma.

Apart from these adverse effects there are always possible complications associated with any surgical procedure such as, but not limited to, infection, nerve damage, and pain which may not be related to the implant.

If an adverse effect occurs, healthcare professionals and/or the patient are obligated by law to inform the manufacturer and patient's home country's National Healthcare authority/ Medicine Agency.

### **Warnings**

- Internal fixation devices aid the surgeon in the augmentation and stabilisation of the cranial bone for osteotomies or augmentation and reconstruction procedures. While these devices are generally successful in attaining these goals, they cannot be expected to replace normal healthy bone or withstand the unsupported stress placed upon the device by full load bearing.
- Internal fixation devices are internal supports, or load sharing devices that protect and construct the defected bone area until normal healing occurs. If the bone healing is delayed or incomplete, the implant can be expected to bend, break, or fail. The size and shape of bones and soft tissue place limitation on the size and strength of implants. Implant size and shape determination and placement must be done carefully, especially in locations of thin overlying tissue such as the corner of the eye.
- It is unlikely for the implant to migrate in the treatment of a bone defect. This is because the implant must always be fixed to the cranial bone with screws due to it being installed on the bone matter.
- Surgical implants are subject to repeated stress in use, which can result in fatigue fracture.
- If the healing of a bone is delayed, unsuccessful or incomplete, the implant may migrate.
- Patient-dependent factors such as each patient's activity level and adherence to loading instructions have an effect on the attachment and/or migration of the implant.
- The surgeon must be thoroughly knowledgeable not only of the medical and surgical aspects of the implant but also the mechanical properties of surgical implants.
- Correct placement of the implant is of great importance. Improper placement may harm surrounding tissues, for example the dura mater. The skin above the implant must be in good condition before surgery and incision above the implant should be avoided.
- Correct handling of the implants is extremely important. The implants should not be shaped or bent so much as to break them. Bending, scratches and notches add to the risk of particle release and implant breakage.
- The manufacturer does not recommend cutting the plate, because this may result in the plate breaking and particles being released from the plate. Dura mater tack-up sutures, for example, must be attached using the premade holes in the implant. Additional holes in the locations desired by the surgeon can be made in the manufacturing phase. If necessary, it is possible to drill additional screw holes at the edges of the plate, but this may lead to small amounts of particles being released from the plate. Any plate modification by cutting or drilling must be carried out in a manner which prevents loose particles from entering the patient's body. Under no circumstances may additional holes be made in areas containing bioglass particles. If additional holes are necessary for dura mater tack-up sutures, the holes must be made at the positions of the reinforcement strips. The manufacturer recommends using the premade holes in the plate, as regards fixation holes and dura mater tack-up sutures. If the need for additional holes presents itself during surgery, they must be implemented in accordance with the above-mentioned instructions.
- Implants can loosen, fracture, migrate, or cause pain. If there is need for implant removal, the implants can be removed after healing. However, this can be challenging as the material might be strongly attached to the bone. Removal might require drilling of the material. The surgeon should weigh the risks versus benefits when deciding whether to remove the implant. Adequate postoperative management should follow implant removal.
- Adequately instruct the patient. Postoperative care is important. The patient's ability and willingness to follow instructions is one of the most important aspects of successful treatment. Patients with senility, mental illness, alcoholism, or drug abuse may be at higher risk of device failure since these patients may ignore instructions and activity restrictions.

- The patient is to be made fully aware and warned that the device does not replace normal healthy bone, and that the device can break, bend or be damaged as a result of stress, activity, load bearing or an accident directed to the device. The patient is to be made aware and warned of general surgical risks, complications, possible adverse effects, and to follow the instructions of the treating physician. The patient is to be advised of the need to come in for regular postoperative follow-up examinations for as long as they are necessary to ensure recovery.

#### **Precautions**

Do not reuse implants. While an implant may appear undamaged, previous stress may have created imperfections that would reduce the service life of the implant. Do not treat patients with implants that have been even momentarily placed in a different patient. Use a GLACE™ Custom Made CMF Implant only on the patient, for whom it is designed. Pay attention to the correct positioning of the plate according to preoperative imaging. GLACE™ Custom Made CMF Implants are MR safe.

#### **Materials**

Glass fibre reinforced polymer structure. Bioactive glass granules, 28–41 w-%, are incorporated into the space between the laminates of the implant. The polymer is dimethacrylate resin polymerised to a high degree of monomer conversion.

#### **Accessories**

- The GLACE™ Custom Made CMF Implant is attached to bone with self-drilling titanium fixation screws with thread diameter of 1.5 mm and appropriate length, for example threaded length of 4–6 mm.
- 3D-model of patient's skull (SIC-1002) and template of the implant (SIC-1001) can be delivered as part of the order, if requested. These can be used as tools to plan the surgery. The template is a non-sterile instrument. The template must never be implanted.
- A GLACE™ Drawing Guide (SIC-1003) can be delivered as part of the order, if requested. This can be used as a tool to plan resection. The drawing guide is delivered STERILE. In case of Drawing Guide has been delivered same time, read its instructions for use before starting the operation. The drawing guide must never be implanted.

#### **Sterility**

The GLACE™ Custom Made CMF Implant is supplied STERILE and must not be re-sterilised. Do not use the implant if the sterile packaging is damaged. GLACE™ Custom Made CMF Implants are sterilised using hydrogen peroxide. The GLACE™ Drawing Guide is supplied STERILE and must not be re-sterilised. Do not use the drawing guide if the sterile packaging is damaged. The GLACE™ Drawing Guide is sterilised using hydrogen peroxide. The template of the implant is delivered NON-STERILE.

#### **Storage**

Do not use if this product has not been stored in accordance with the following storage conditions:

- Temperature: 15–28 °C
- Humidity: 0–80% relative humidity, non-condensing

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Skulle Implants Oy